



Capture critical biomarkers from challenging cancer samples




Low input, high-/low-quality, FFPE RNA-seq solutions




High-quality libraries from a range of samples types and quality (RIN 3-10)

- Including fresh, fresh frozen, and highly degraded FFPE samples
- Compatible with even < 1 ng RNA



Capture coding and non-coding transcripts with high sensitivity & specificity

- Detect mRNA, non-PolyA RNA, long non-coding RNA (lncRNA) species
- Prepare stranded libraries



All-in-one kit for library prep and ribosomal RNA depletion

- Streamlined, single day protocol ~7h



Improve data accuracy with unique molecular identifiers (UMIs) and unique dual indexes (UDIs)



Learn more at
takarabio.com/totalrna-seqkits

that's
GOOD
science!

RNA biomarker analysis

High-quality data from even highly degraded FFPE samples across a wide range of inputs




Figure 1. Evaluation of input RNA integrity and NGS library profiles for FFPE samples. Panels A–D. Bioanalyzer traces of RNA inputs obtained from four different samples: one healthy breast tissue sample, and three lung tissue samples obtained from cancer patients. RNA profiles and DV200 values indicate the integrity of each sample. Panels A'–D'. Bioanalyzer traces of sequencing libraries generated from the corresponding RNA inputs. The profile data suggests that high-quality libraries are produced regardless of input RNA integrity.




Figure 2. Sequencing metrics for FFPE samples show comparable numbers of transcripts identified with fragments per kilobase per million reads mapped (FPKM) >1 , and a high degree of correlation across input amounts.




Figure 3. Comparison of transcript expression measurements shows highly reproducible transcript expression across input amounts for samples of varying integrity including two highly degraded samples with DV200 $\sim 30\%$ (C and D). Plots A – D show data from samples with corresponding labels in Figure 1.

Choose the right kit

	SMARTer® Stranded Total RNA-Seq Kit v3 - Pico Input Mammalian	SMARTer® Stranded Total RNA-Seq Kit v2 - Pico Input Mammalian	SMART-Seq® Stranded Kit
Input	250 pg–10 ng, total RNA, 10–1,000 cells	250 pg–10 ng, total RNA	10 pg–10 ng, total RNA, 1–1,000 cells
FFPE compatible	Yes	Yes	No
RNA integrity	RIN 3–10; DV ₂₀₀ $>25\%$ (including FFPE)	RIN 3–10; DV ₂₀₀ $>25\%$ (including FFPE)	RIN 3–10
rRNA removal	H/M/R* rRNA	H/M/R rRNA	H/M/R rRNA
UDIs	384 UDIs	384 UDIs	384 UDIs
UMIs	Yes	No	No
Workflow time	~ 7.5 h	~ 6 h	~ 7 h

*H/M/R = human, mouse, rat



Learn more at
takarabio.com/totalrna-seqkits

that's
GOOD
science!