

# CERTIFICATE OF ANALYSIS - TECHNICAL DATA SHEET

**Product name** C5a, Mouse, Recombinant

Catalog number HC1101

Lot number - Expiry date -

**Volume** Reconstitute with 0.5 ml distilled/de-ionized water **Activity** N.A.

Formulation Lyophilized in 0.005% Tween 80+2mM reduced Amount

L-glutathion+0.2mM oxidized

L-glutathion+0.1M Tris-Cl buffer, pH 8.0

Host SpeciesMouse, E.coli-derivedConcentrationN.A.

Endotoxin N.A. Purification N.A.

Storage 4°C Purity N.A.

## **Application notes**

For dilutions use protein stabilized phosphate buffered saline, pH 7.4.

#### **General Information**

### Description

Complement fragment C5a is a 74-residu glycopolypeptide which is generated by proteolytic cleavage of the complement factor C5 in the course of complement activation. C5a is a potent chemoattractant and anaphylatoxin that acts on all classes of leukocytes and on many other cell types including endothelial, smooth muscle, kidney, liver, and neural cells. In addition to its proinflammatory effects, C5a has been shown to protect cells against toxic insult and to stimulate proliferation in neurons and hepatocytes, suggesting a wider role for C5a in homeostasis. C5a is rapidly desarginated by serum carboxypeptidase N to the less potent derivate C5a desArg, the first stage in deactivation of anaphylatoxin activity. The C5a desArg form has a different spectrum of bioactivity to intact C5a. Recombinant mouse C5a is His-Tagged and has the following amino acid sequence:

>50µg

MRGSHHHHHHGSDYDIPTTENLYFQGGSNLHLLRQKIEEQAAKYKHSVPKKCCYDGARVNFYETCEERVARVTIGPL CIAFNECCTIANKIRKESPHKPVQLGR. The molecular weight of the protein is 12 kg/mol.

## Aliases

Pam3Cys-Ser-(Lys)4 x 3 HCl, Pam3Cys-SK4 x 3 HCl, P3CSK4 x 3 HCl

### References

- Bosmann, M et al; The interaction between C5a and both C5aR and C5L2 receptors is required for production of G-CSF during acute inflammation. Eur J Immunol 2013, 43:1907
- Ma, N et al; C5a Regulates IL-12+DC Migration to Induce Pathogenic Th1 and Th17 Cells in Sepsis. PLoSONE 2013, 8:e69779
- Andersson, C et al; Rapid-onset clinical and mechanistic effects of anti-C5aR treatment in the mouse collageninduced arthritis model. Clin Exp Immunol 2014, 177:219

### Storage&stability

Lyophilized product should be stored at 4°C. Store stock solution after reconstitution in aliquots at -20°C. Repeated freeze and thaw cycles cause loss of activity. Under recommended storage conditions, product is stable for at least one year.

#### **Precautions**

Caution: vial is under vacuum. For research use only. Not for use in or on humans or animals or for diagnostics. It is the responsibility of the user to comply with all local/state and federal rules in the use of this product. Hycult Biotech is not responsible for any patent infringements that might result from the use or derivation of this product.

We hereby certify that the above-stated information is correct and that this product has been successfully tested by the Quality Control Department. This product was released for sale according to the existing specifications. This document has been produced electronically and is valid without a signature.

Approved by Manager of QC Robbert Zwinkels

Date 23/03/2018

Version: 02-2018

Do you have any questions or comments regarding this product? Please contact us via <a href="mailto:support@hycultbiotech.com">support@hycultbiotech.com</a>.